STA[®] - CaCl₂ 0,025 M
Solution de CaCl₂ 0,025 M pour les tests de coagulation

- Coffret contenant 24 flacons de 15 ml

REF: 00367)  

2940 07 - Mai 2014 Français 1

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REF: 00367)  

May 2014 English

1/ INTENDED USE
STA[®] - CaCl₂ 0,025 M is an aqueous solution of calcium chloride intended for the coagling tests such as the activated partial thromboplastin time (APTT) or for the assays of intrinsic pathway factors.

2/ KIT CONTENTS
An Asepy Value insert with a barcode is provided in the box. This barcode contains the following information: lot number, kit code number, reagent code number and expiration date.
STA[®] - CaCl₂ 0,025 M : 0,025 M CaCl₂ solution.
This solution contains sodium azide (< 1 g/l) as a preservative. Solutions containing sodium azide should be discarded with care to prevent the formation of explosive metallic azides. If waste materials are dumped into sinks, use copious quantities of water to flush plumbing thoroughly.

3/ CAUTION
Store at 2-25 °C. For *in vitro* diagnostic use only. This solution is to be used only by certified medical laboratory personnel authorized by the laboratory. When this STA[®] - CaCl₂ 0,025 M is used with analyzers of the STA[®] line, read the Reference Manual of the analyzer model carefully before starting. Exercise great care in the handling of this solution and of patient samples. The disposal of waste materials must be carried out according to current local regulations.
In the USA, wherever appropriate, observe CLIA-88 requirements.

4/ SOLUTION PREPARATION
If the solution is refrigerated, allow it to stand at room temperature (18-25 °C) for 30 minutes, before use.
Do not install either an STA[®] - Reducor or a perforated cap on the solution bottle if it is to be used on analyzers of the STA[®] line.
5/ SOLUTION STORAGE
The solution in intact bottles is stable until the expiration date indicated on the box label, when stored at 2-25 °C.
After opening, it remains stable for:
– 24 hours at 37 °C
– 3 days on STA-R[®] and STA Compact[®]
– 6 days on STA Satellite[®]
In case of partial use, the remaining solution stored at 2-8 °C in its original bottle with the cap on top, is stable, when free of any contamination, until the expiration date indicated on the bottle.

6/ REAGENTS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Reagents for performing coagulation tests.
- Analyzer such as STA-R[®], STA Compact[®] or STA Satellite[®].
- Common clinical laboratory equipment and materials.

7/ PROCEDURE

- When performing tests with semi-automated analyzers, refer to the test reagents instructions for use.
- With analyzers of the STA[®] line, after opening of the solution bottle, load it into the instrument according to the recommendations of the Reference Manual of the analyzer model. The vial position in the instrument is the following:
– on STA-R[®] model, place the bottle in the R2 area of the product drawer
– on STA Compact[®] model, place the bottle in the product drawer
– on STA Satellite[®] model, place the bottle in the product carousel.
The 0,025 M CaCl₂ solution is automatically used by the instrument (see the Reference Manual).

Significant changes are indicated by dotted lines in the margin.

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 Information and/or picture contained in this document are protected by copyrights and other intellectual property rights. © 2014, Diagnostica Stago, all rights reserved. No part of this document may be reproduced without the prior written permission of Diagnostica Stago.

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DIAGNOSTICA STAGO S.A.S.