

CAPILLARYS
TERA **3**


Efficiency in lab automation solution

Multi-assay analyzer designed to fit your lab routine

Automated
management




High performance
and autonomy



Simplified reagent
management



Optimized sample
handling



High modularity
to meet your needs

CAPILLARYS **3**

Integrability for
efficient workflow



sebia
The new language of life

Reduction in operator intervention *for high quality results delivered on time*

High capacity for on time result delivery

- Excellent resolution for high analytical performance
- High throughput

CAPILLARYS 3 Test Menu

Hb A1c	63 tests/h
Serum Protein Electrophoresis	99 tests/h
Urine Protein Electrophoresis	97 tests/h
Serum Immunotyping	11 tests/h
Urine Immunotyping	11 tests/h
CDT/CDT _{IFCC}	73 tests/h
Hemoglobin Electrophoresis	65 tests/h

For up to date information, refer to your local representative.



Simplified reagent management

- High on-board reagent and waste capacity
- On-the-fly reagent change
- RFID tags for reagents traceability



Automated instrument management

- LCD touch screen to simplify the instrument operations
- Automated and programmable maintenances
- Simplified Quality Control management



Optimized sample handling

- High sample loading capacity of various tube types
- Capped and uncapped tubes
- Sample mixing capability



Specifications


- Length: 90 cm • Height: 54 cm
- Depth: 67 cm • Weight: 75 kg

Ordering *information*

Instrument

- CAPILLARYS 3 TERA1246

For further information, please find your Sebia local representative to contact on our Website
www.sebia.com/sebia/worldwide-presence/.



Scan here
for more
information on
CAPILLARYS 3 TERA



Sebia | RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008 Evry cedex | France | Tel.: +33 1 69 89 80 80 | E-mail: sebia@sebia.com | www.sebia.com

©2025  Sebia - RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008 Evry cedex | France - Ref.: C3TFIEN - 01/2023
Non-contractual illustration - Sebia reserves the right to modify at any time and without notice the information contained in this document intended for healthcare professionals. This document contains information intended for wide distribution and may therefore contain product details or information that is not available or valid in your country. Carefully read the instructions in the reagent package inserts and instrument manuals.